





Learning Objectives



- Terminal Learning Objective



- **Old**

- 1-Eye Patch Dressing
- 1-Battle Dressing with Tie Ends
- 1-Bottle of Iodine Water Purification Tablets
- 1-Bottle of Wound Disinfectant
- 1-Lip Balm
- 2-Cravat's
- 10-Adhesive Bandages



- **New**


- 1-Packet of Adsorbent Hemostatic Agent (QuikClot)
- 1-Package of Burn Gel (Water-Jel)
- 1-Bottle of Iodine Water Purification Tablets
- 1-Bottle of Wound Disinfectant
- 1-Tourniquet
- 2-Pressure Dressings
- 2-Bulky Gauze Rolls
- 10-Adhesive Bandages
- 5-Large Adhesive Bandages
- 2-Triangle Bandages



Components and Characteristics



- 1-Packet Adsorbent Hemostatic Agent
 - Granular mineral powder
 - Sterile
 - Stops high volume venous and arterial bleeding through rapid absorption of liquid
 - Clot remains in the wound until removed by medical personnel
 - Food and Drug Administration (FDA) approved





Components and Characteristics (Cont)



- 1-Package 4"x16" Sterile Burn Dressing (Water-Jel)
 - Extinguishes white phosphorous burns
 - Cools burn area
 - Reduces chance for hypothermia
 - Prevents burn from further progression
 - Protects against infection and eases pain
 - FDA approved
 - Do not use for other than intended use





Components and Characteristics (Cont)



- 1-Bottle of Iodine Water Purification Tablets
 - To purify non potable water for consumption
 - After handling the tablets, wash hands thoroughly. If eyes come in contact with the tablets or residue, it could cause temporary blindness





Components and Characteristics (Cont)



- 1-Bottle of Wound Disinfectant
 - Betadine solution
 - Cleansing of Minor Wounds





Components and Characteristics (Cont)



- 1-Tourniquet
 - One Handed Use
 - Use on arterial and venous blood
 - For Moderate to Severe Bleeding only
 - For use on Extremities only





Components and Characteristics (Cont)



2-Pressure Dressings

- Vacuum Packed
- 8" x 10" Absorbent Pad
- 4" Wide elastic wrap with Velcro strip on both ends
- Coated Steel tension hook for individual application





Components and Characteristics (Cont)



- 2-Bulky Gauze Rolls
 - Gauze roll
 - Absorbent Material



- 10-Adhesive Bandages





Components and Characteristics (Cont)



- 2 Non-Sterile Triangular Bandages
- 40" x 40" x 56"





Questions?